

# Development of Mask Contamination/Inspection System by coherent EUV light

Sangsul Lee<sup>1,2</sup>, Chang Young Jeong<sup>1</sup>, Seungyu Rah<sup>3</sup>, Donggun Lee<sup>4</sup>,  
Seong-Su Kim<sup>4</sup>, Han-Ku Cho<sup>4</sup> and Jinho Ahn<sup>1\*</sup>

<sup>1</sup> Department of Materials Science and Engineering, Hanyang University

<sup>2</sup> Information Display Research Institute, Hanyang University

<sup>3</sup> Pohang Accelerate Laboratory (PAL)


<sup>4</sup> Photomask Team, Memory Division, Semiconductor Business, Samsung Electronics Co., LTD



- EUVL Test Bed – Infrastructure for development of EUVL
  - Introduction of EUVL beamline @ Pohang Light Source.
- Mask Contamination/Inspection System
  - What is it's concept?
  - Optical design for mask contamination/inspection system
  - Test results
    - : In-situ mask contamination analysis
    - : mask inspection - Light source alignment & Noise reduction
- Summary



## Development of the Core EUVL technology for 22 nm node



## EUV Focus Areas 2003-2008




2003	2004	2005	2006	2007	2008
1. Source power and lifetime including condenser optics lifetime	1. Availability of defect free mask	1. Resist resolution, sensitivity & LER met simultaneously	1. Reliable high power source & collector module	1. Reliable high power source & collector module	1. Long-term source operation with 100 W at IF and 5MJ/day
2. Availability of defect free mask	2. Lifetime of source components & collector optics	2. Collector lifetime	2. Resist resolution, sensitivity & LER met simultaneously	2. Resist resolution, sensitivity & LER met simultaneously	2. Defect free masks through lifecycle & inspection/review infrastructure
3. Reticle protection during storage, handling and use	3. Resist resolution, sensitivity & LER met simultaneously	3. Availability of defect free mask	3. Availability of defect free mask	3. Availability of defect free mask	3. Resist resolution, sensitivity & LER met simultaneously
4. Projection and illuminator optics lifetime	• Reticle protection during storage, handling and use	4. Source power	4. Reticle protection during storage, handling and use	4. Reticle protection during storage, handling and use	• Reticle protection during storage, handling and use
5. Resist resolution, sensitivity and LER	• Source power	• Reticle protection during storage, handling and use	5. Projection and illuminator optics quality & lifetime	5. Projection and illuminator optics quality & lifetime	• Projection / illuminator optics and mask lifetime
6. Optics quality for 32-nm half-pitch node	• Projection and illuminator optics lifetime	• Projection and illuminator optics quality & lifetime			




Source: 2008 EUVL Symposium



# EUVL Test Bed – Infrastructure for EUVL



Recent performance.  
ICS/CSM to investigate mask life time issue related with C contamination




# EUVL Test Bed – Recent progress

Unit	Resolution	Photo-speed (@50nm, L/S)	LER (@50nm)
07' 4Q	< HP 32 nm	20.0 mJ	5.9 nm
08' 2Q	<b>HP 28 nm</b>	<b>6.2 mJ</b>	<b>5.3 nm</b>


## LBNL test results

**Exposure : LBNL MET**  
**Illumination : Ann (0.3/0.55)**  
**Inspection : L/S 1:1**  
**Film Thickness : 50nm**

## Dongjin DHE-1158



2008.05. LBNL test result



- Scan Speed : # 1
- Baseline Noise :  $7 \times 10^{-15}$
- Scan Time : 10 min
- Beam Current : 441 nA


Resist outgassing test result  
 -> Dongjin EUV resist evaluation




# EUVL Test Bed – new challenge (mask inspection)

1. EUV mask analysis (CD, H-V bias)
2. In-situ mask contamination & analysis (reflectance, carbon thickness)

**In-situ** measurement of mask contamination effect on the imaging performance





## EUV CSM


- EUV coherent scattering microscopy
- Mask image reconstruction using EUV diffraction beam
- Critical dimension (CD) measurement

## In-situ contamination system

- Exposure acceleration system
- ML reflectance & carbon thin film thickness measurement(In-situ spectroscopic ellipsometer)
- Exposure environment optimization for lower contamination





# Mask Contamination/Inspection System




X-ray CCD camera  
: Pixel - 2048 X 2048 imaging array  
Pixel size - 13.5 X 13.5  $\mu\text{m}^2$   
Imaging area - 27.6 X 27.6  $\text{mm}^2$   
Vacuum compatible


microscope  
: 85mm beam  
distance  
X 5000 zoom



# Optical design for Mask Contamination/Inspection System




# Optical design for Mask Contamination/Inspection System



	Optics Name	Size	Radius of Curvature(mm)	Beam incidence angle to optics (deg.)
1	Pinhole	5 (um)	inf.	0
2	Spherical mirror	1"	165	3
3	Plane mirror	1/2"	inf.	39
4	Mask	6"	inf.	6
5	CCD center	1"	inf.	0
6	Plane mirror for imaging	1/2"	inf.	26
7	Spherical mirror for imaging	1/2"	88	16
8	Mirror for acceration	1/2"	inf.	42
9	Mask	6"	inf.	6
10	Reference mirror	1"	inf.	6
11	Photodiode-2nd(sqar)	10 mm	inf.	0
12	Photodiode-1st(quad)	10 mm	inf.	0




# Mask Contamination Analysis




1. Reference ML mirror
2. EUV mask acceleration test (XX hours) : Change of EUV reflectance  
Carbon thin film thickness


## \* EUV acceleration, 2hr



**EUV reflectometer**




Blank mask (A) : 1.85% Reflectivity Drop



Blank mask (B) : 2.92% Reflectivity Drop



# Mask Contamination Analysis - ellipsometry





Carbon contamination image  
(Measurement of In-situ microscope)

: Si capping mirror

EUV expose : 2hr

Thickness : ~420 Å

$n : 1.6 \quad k : 0.0 \quad @ \lambda = 633\text{nm}$




50 $\mu\text{m}$  ~ 1 mm beam diameter

In-situ measurement  
: carbon thickness

65° incident angle




# Development of Coherent Scattering Microscopy





## EUV Coherent Scattering Microscopy

Mask Pattern (SEM Image)



100-nm (4X)





Field Spectrum Measurements


Source : EUV symposium 2008  
Donggun Lee, et al




# Beam alignment & noise reduction



remove  
noise  
source



back  
ground  
subtraction



We set up Mask Contamination/Inspection system at PAL.

Our system will help to quantize and correlate the degree of contamination and the degradation of mask performance.



# Acknowledgements

The authors would like to thank the team of Korean EUVL society.

This work was supported in part by the EUVL R&D research fund of MOCIE and internal research support program of New Growth Engine Semiconductor Center.

